
Invacare® **Etude**

User's manual (DK - S - GB - NL - D - F - ES - P)

Invacare® EC-Høng AIS
Ident. nr.: I423886
Version 01 06.2002

Invacare® **Etude**

Invacare® Etude

BRUGSANVISNING (DK)	3 - 17
BRUKSANVISNING (S)	18 - 32
USER'S MANUAL (GB)	33 - 47
HANDLEIDING VOOR GEBRUIK (NL)	48 - 62
BEDIENUNGSANLEITUNG (D)	63 - 77
MANUEL D'UTILISATION (F)	78 - 92
INSTRUCCIONES PARA EL USO (ES)	93 - 107
MANUAL DO UTILIZADOR (P)	108 - 122

Indholdsfortegnelse

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Table of Content

1. General information	34
2. Unpacking your product	36
3. Assembly of Etude	37
4. Fitting the accessories	39
5. Operating the bed	39
6. Emergency lowering of the backrest and/or thigh section ..	40
7. Operating the accessories	40
8. Dismantling the Etude	41
9. Using the transporter brackets	42
10. Order numbers for accessories	42
11. Cleaning and disinfection	42
12. Maintenance and check-ups	42
13. Maintenance chart	43
14. Waste disposal	43
15. Trouble-shooting with the electrical system	44
16. Wiring	45
17. Technical specifications	47
18. Weight	47

Congratulations

You have chosen an Invacare® Etude nursing bed. Etude is the detachable nursing bed from Invacare® EC-Høng designed especially for home nursing. The Etude nursing bed provides ease of handling combined with good functionality.

Invacare® EC-Høng is certified according to DS/EN ISO 9001/EN 46001 which ensures that our customers are always supplied with products of uniform quality.

Throughout the entire production process our materials/products are quality controlled by the operators. Moreover, a final test is carried out prior to packaging and shipping.

The operator carrying out the final test, which comprises of checking all moveable parts, motors and castors, issues a QA number on the product thereby confirming its quality.

QA XXX

If the product does not correspond to the quality demands of Invacare® EC-Høng, it will be discarded.

If, contrary to our expectations, a problem should arise in connection with the delivered product, please contact your Invacare® supplier.

Invacare® reserves the right to alter product specifications without prior notice.

Invacare® accepts no liability for any use, change or assembly of the product other than as stated in this User's Manual.

Accessories not mentioned in this manual must not be used.

I. General information

- The Etude bed is marked with the  mark in accordance with directive 93/42/EEC concerning Medical Apparatus.
- The Etude bed has been tested and approved according to EN 1970:2000.
- The motors and control of the Etude bed have been approved according to EN 60601:1996-03.
- The Etude bed has been approved and marked with the TÜV mark.
- Beds provided with tilt function cannot be approved according to EN 1970:2000.
- The Etude bed has undergone a risk analysis according to EN 1441.
- The control and motors are protected according to IP 54. The hand control box is protected according to IP 66.
- Max. weight: Etude Basic: 180 kg - Etude Classic: 220 kg.
- Max. patient weight: 145 kg (Etude Basic) - 185 kg (Etude Classic) provided that the weight of the mattress and the accessories do not exceed 35 kg.
- Standard size of mattress for the Etude Basic is: 90x200x12, Etude Classic: 90x200x10.
- If the patient height exceeds 2 meters it is recommended to use a mattress support extension.
- The Etude bed has loose cables, which have to be exchanged when damaged.
- The bed is not intended for children under 12 years and psychiatric patients.
- Remove the plug from the mains voltage before moving the bed. The cable must be kept clear off the floor and the castors during transportation.
- The adjustment area of the mattress support is: Basic: 40-80 cm or 33-73 cm, Classic: 40-80 cm (30-70 cm/Low height version).
- Only accessories described in this manual may be mounted on the Etude bed.
- The angle between the lower leg section and horizontal is adjustable from 0° to 15°.

Note! The total weight (mattress, accessories and patient weight) must not exceed 180 kg (Basic), 220 kg (Classic).

The bed is not intended for use together with MEDICAL ELECTRICAL EQUIPMENT.

It can be dangerous to roll over the mains cable.


Do not bring mains cable into moving parts.

Disconnect the plug from the mains before moving the bed.

The cables must be mounted in such a way that they are kept clear off the floor and do not block the castors.

We recommend to mount the mains cable on the hook for this purpose, see picture below.

The bed's castors must be locked during nursing a patient in bed or when using any positioning function. The bed must not be moved when occupied.

Adjust the mattress support sections to horizontal position and lower the bed to its lowest position before moving the bed.

Hold the top of the bed end with both hands while the bed is pushed/pulled.
Any obstacle must not be higher than 10 mm.

If the functions of the bed change, check the bed according to the maintenance chart, chapter 13.

All service and maintenance work described in chapters 11, 12 and 13 must only be performed by personnel who have been instructed and trained by Invacare[®].

The Etude bed must be stored in a room with a humidity equivalent to normal indoor conditions.

The bed must be stored in a frost-free environment.

2. Unpacking your product

For a complete bed the following parts must be provided:

(Illustration 1, here shown as a 2-sectioned bed)

1. Etude Basic

A) 1 mattress support

If 2-sectioned = ETUDE.CXXXX.XX (see complete list for specific number)

If 4-sectioned = ETUDE.EXXXX.XX (see complete list for specific number)

2 mattress handles, hand control and motors are mounted on the mattress support.

B) 1 pair of bed ends = ETUDE.000CI.XX

Optional:

- Pair of siderails (steel)
- Lifting pole
- Support handle(s)
- Mattress support extension (10 or 20 cm)
- Transport fitting

2. Etude Classic

A) 1 mattress support

If 2-sectioned = ETUDE.CXXXX.XX (see complete list for specific number)

If 4-sectioned = ETUDE.EXXXX.XX (see complete list for specific number)

2 mattress handles, hand control and motors are mounted on the mattress support.

B) 2 individual bed ends

If standard = ETUDE.000EI.XX

If low = ETUDE.000FI.XX

Optional:

- Pair of siderails (wood)
- Lifting pole
- Support handle(s)
- Mattress support extension (20 cm)
- Transport fitting

Order numbers for Etude mattress support and bed ends:

ETUDE.CE100.XX 2-SECT 90X200 STD. MAINS CABLE, TILT, NOT TÜV APPROVED

ETUDE.CE300.XX 2-SECT 90X200 STD. MAINS CABLE, EXCL.TILT, TÜV APPROVED

ETUDE.CU100.XX 2-SECT 90X200 UK MAINS CABLE, TILT, NOT TÜV APPROVED

ETUDE.CU300.XX 2-SECT 90X200 UK MAINS CABLE, EXCL.TILT, TÜV APPROVED

ETUDE.CA100.XX 2-SECT 90X200 AU/NZ MAINS CABLE, TILT, NOT TÜV APPROVED

ETUDE.CA300.XX 2-SECT 90X200 AU/NZ MAINS CABLE, EXCL.TILT, TÜV APPROVED

ETUDE.EE100.XX 4-SECT 90X200 STD. MAINS CABLE, TILT, NOT TÜV APPROVED

ETUDE.EE300.XX 4-SECT 90X200 STD. MAINS CABLE, EXCL.TILT, TÜV APPROVED

ETUDE.EU100.XX 4-SECT 90X200 UK MAINS CABLE, TILT, NOT TÜV APPROVED

ETUDE.EU300.XX 4-SECT 90X200 UK MAINS CABLE, EXCL.TILT, TÜV APPROVED

ETUDE.EA100.XX 4-SECT 90X200 AU/NZ MAINS CABLE, TILT, NOT TÜV APPROVED

ETUDE.EA300.XX 4-SECT 90X200 AU/NZ MAINS CABLE, EXCL.TILT, TÜV APPROVED

ETUDE.EJ100.XX 4-SECT 90X200 JPN. MAINS CABLE/100V, TILT, NOT TÜV APPROVED

ETUDE.EJ300.XX 4-SECT 90X200 JPN. MAINS CABLE/100V, EXCL.TILT, TÜV APPROVED

All the mattress supports can be combined with these 3 available types of bed ends:

BED ENDS BASIC (Ordered and delivered as a pair of bed ends)

ETUDE.000CI.XX | **PAIR** BASIC BED ENDS, 40-80 CM (33-73)

BED ENDS CLASSIC (Ordered and delivered as 1 pcs. of bed end)

ETUDE.000EI.XX | **PCS.** CLASSIC BED END, 40-80 CM

ETUDE.000FI.XX | **PCS.** CLASSIC BED END **LOW HEIGHT**, 30-70 CM

XX is colour code.

For standard light grey (RAL 7032) **.21** is to be used for bed ends and **.20** for mattress support.

Example:

ETUDE.000EI.21 (1 Etude Classic bed end in light grey)

ETUDE.EE100.20 (1 Etude 4-sect. mattress support in light grey)

3. Assembly of Etude

(Illustration 2a, 2b+c, 2d)

a) Mattress Support

- Place the two inserts in the head end of the mattress support. The inserts must be mounted in such a way that one extends further out of the side tube than the other.
- Loosely screw in the two thumb screws.
- Push the foot end of the mattress support onto the two inserts.
- Tighten with two thumb screws.
- Retighten the two thumb screws at the head end of the top frame.
- Turn the two mattress handles of the backrest up.
- Turn the handle of the leg section up.

b) Etude Basic

- Turn the locking rings*) on the bed ends to the position "open".
- Latch the mattress support to the bed ends and press firmly into position.
- Turn the locking rings to the position "locked".

c) Etude Classic

- Release the quick locks*) on the bed ends.
- Latch the mattress support to the bed ends and press firmly into position.
- Turn the quick locks to the position "locked".

*) The locking rings/quick locks ensure that the mattress support is not accidentally lifted clear of the bed ends.

d) Control Box

The control box is latched to the backrest motor.

The control box is provided with a label with symbols showing where to connect the motors' plugs.

- 1: Backrest motor
- 2: Thigh section's motor
- 3: Bed end
- 4: Bed end
- HB: Hand control

Wiring (Also see the pictures in section 16)

Attention!

In order to prevent the cables from being torn apart when activating the motors it is critical that you follow the instructions below:

- 1) Connect the cables of the bed end motors (head and foot end) directly to the control box.
- 2) Only 4-sectioned beds: The plug of the thigh motor must be directed through the opening of the supporter for the backrest motor and then be connected to the control box.
- 3) Connect the mains cable to the 230V plug.
- 4) Run the motors of the bed ends to their top position.
- 5) Place the motor cable of the foot end on the hook at the foot end.
- 6) Run the backrest to its top position.
- 7) Place the head end motor cables on the 4 hooks at the head end.
- 8) Fasten the safety comb to the control box.

When working with the foot end, note that there is a risk of getting squeezed during transport since the adjustable leg section is not firmly locked. Remember to remove the plug from the mains socket when moving the bed.



4. Fitting the accessories

(Illustration 3)

Fitting the siderail, model 5519/5525/5526/5529

The siderail must be mounted with the release system at the head of the bed.


The fork links on the siderail must be mounted according to the instruction on the siderail.

Tighten the siderail with 2 thumb screws.

Mounting of BRITT II/BRITT III (Classic bed ends)

1. Raise the bed about 1/3.
2. Dismount the stop-knobs using the thumb screws (C).
3. Remove protection tape at the locking pins at both ends of the siderail.
4. Fit one end at a time. Lift the lower wooden bar and guide the metal bow to the bed end rail (A).
5. Press the locking pin of the siderail with a finger (B).
6. Guide the siderail to the bed end and pull up the wooden siderail bar, until the locking pin locks with an audible click.
7. Mount the stop-knobs using the thumb screws (C).

There is a risk entrapment on assembly and operation of the siderails.

Fitting the lifting pole

Remove the plug from the lifting pole tube at the head end of the bed where the lifting pole is to be placed. Insert the lifting pole into the lifting pole tube and fix with the thumb screw.

5. Operating the bed

Electric operation

Manual operation

Adjustment of backrest (Illustration 4) Raising or lowering the backrest to the required angle: Use the button with the symbol shown on the right.		
Adjustment of the leg section (Illustration 5) Use the handle when raising or lowering the leg section. Raising or lowering the leg section to the required angle: UP: Lift the leg section. DOWN: Lift the leg section right up and then lower it.		
Adjustment of the thigh section (only 4-sectioned beds) (Illustration 6) Use the button with the symbol shown on the right.		
Adjustment of height of mattress section (Illustration 7) Use the button with the symbol shown on the right.		

Always leave the bed in the lowest position. Otherwise there is a risk of entrapment due to accidental lowering of the mattress support. A person under the bed can be seriously injured during height adjustments.

There is a risk of entrapment between the thigh section and the mattress support when the thigh section is lowered.

There is a risk of entrapment between the backrest and the cross rail on the head end of the mattress support. There is also a risk of getting squeezed between the leg section and the cross rail on the foot end of the mattress support when the backrest/leg section is lowered.

If the bed is equipped with tilt function, please note that the bed must not be used without the supervision of the nursing staff.

Brakes

(Illustration 8)

Each of the bed's 4 wheels is provided with brakes for both lengthwise and crosswise locking. The brake is foot-operated.

6. Emergency lowering of the backrest and/or thigh section

(Illustration 9)

An emergency release of the mattress support may be necessary in case of e.g. a power or motor failure. An emergency release of the height adjustment is NOT possible!

- Remove the plug from the mains before emergency release of the mattress support.
- In an emergency the mattress sections are released by pulling out the cotter pin from the motor.

A minimum of 2 persons is required to release a mattress section. Both persons hold the mattress section in locked position. One of them pulls out the cotter pin. Both slowly lower the mattress section until it is completely down.

7. Operating the accessories

Operation of siderail, model 5519/5525/5526/5529

(Illustration 10)

UP: Pull the top tube on the siderail towards the end where the release system is placed.



DOWN: Press the release button down while pulling the top tube away from the release system.

Operation of siderail, model BRITT II/III

The bed end may be equipped to lock the siderail at half height as well as in the top position. Normally the siderail can only be locked in top position.

UP: Pull up the top wooden siderail bar, until the locking pin locks with an audible click.

DOWN: Lift the top wooden siderail bar and press the two locking rings together. Lower the siderail.

There is a risk of entrapment when operating the siderails.

Adjusting the height of the lifting pole handle

Loosen the cord as shown in fig. A. The lifting handle can now be adjusted to the desired height. Press the cord together as shown in fig. B and check that the cord is locked in the cord lock by pulling the handle.

Position the lifting pole in such a way that the handle extends inwards across the bed. If the lifting pole is used while the handle has been turned away from the bed - the bed can tip when the handle is used.


fig. A Fitting the fig. B

mattress support extension

Disassemble the bed and remove the inserts. Mount the mattress support extension on either head or foot end and reassemble the bed. (See illustration 2).

8. Dismantling the Etude

- Dismantle siderails and lifting pole.
- Bring the bed to its lowest position and adjust all mattress support sections to horizontal position.
- Disconnect the 230 V supply. Roll the cable onto the hook of the head end of the bed.
- Disconnect all the motor plugs from the motors (do not disconnect the backrest motor cable - cable I).
- Dismantle the safety comb from the control by means of a tool, e.g. a screwdriver.
- Pull out the cables of the bed end motors. Only 4-sectioned: Also pull out the cables of the leg section motors.
- Separate the mattress support from the bed ends.
- Divide the head and foot end of the mattress support.

9. Using the transporter brackets

See instructions for use sent separately in the transporter kit boxes.

10. Order numbers for accessories

PLEASE NOTE ONLY ETUDE BASIC IS APPLICABLE TO THE UK.

General (Etude Basic + Etude Classic)

Lifting pole	50.57600.20
Mattress handle	215989
Transport fitting	1423920-7032
Support handle 25x30 cm (1 pcs.)	021964.20
Support handle 25x80 cm (1 pcs.)	1417510-7032
Support handle 40x30 cm (1 pcs.)	1417511-7032
Support handle 40x50 cm (1 pcs.)	021963.20
Support handle 40x95 cm (1 pcs.)	1417512-7032

Etude Basic

Siderail (folding) 1 pcs. right and 1 pcs. left	50.55190.20/L+R
Siderail (folding) 1 pcs. right and 1 pcs. left (25 mm tube)	50.55250.20/L+R
Siderail (folding) 1 pcs. right and 1 pcs. left (with plastic insert)	50.55260.20/L+R
Extended siderail (folding) 1 pcs. right and 1 pcs. left	50.55290.20/L+R
Mattress support extension - 10 cm	1423993-7032
Mattress support extension - 20 cm	1423994-7032

Etude Classic (Not applicable to the UK)

Wooden siderail Britt III (pair)	020434.01
Extended wooden siderail Britt II (pair) +20 cm	019376.01
Mattress support extension - 20 cm	1423994-7032

11. Cleaning and disinfection

Remove the mains plug from the socket before cleaning.

The Etude bed does not tolerate cleaning in an automatic washing plant or using water-jet based cleaning equipment.

The bed should be washed down using a sponge, cloth or brush. Use ordinary household cleaning agents.

When disinfecting only use spray or wiping methods. Only officially recognized disinfectants must be used.

Dry the bed after cleaning.

Never use acids, alkalines or solvents such as acetone or cellulose thinner.

The hand control box, motors and control may be washed with brush and water, but not with pressurized water.

If the backrest is raised be careful not to lower it unintentionally, as there is a risk of entrapment between the backrest and the top frame.

12. Maintenance and check-ups

Service and maintenance of the Etude bed must only be performed by personnel who have received the necessary instruction or training.

Thumbscrews must be checked for tightness and the mechanism of the siderail locking and motion must be checked after 3 months of use.

With normal everyday operation, the instructions of the maintenance chart, chapter 13, must be followed once a year.

Check the bed and particularly the cables, plugs and the condition of the siderails after each use.

13. Maintenance chart

Service and maintenance of the Etude must only be performed by personnel who have received the necessary instruction or training.					
Bed / Id. no.:					
Date:	Initials:				
Visual inspection of all parts of the bed. (Plastic deformation and/or wear and tear of welded joints).					
Control of all centres of rotation (Motors + mattress support parts).					
All motors running without failures. (With regular speed and at low noise).					
Check the rastofix fittings and function.					
Visual inspection of all cabinets for damages.					
Check that mains cable and plug are intact. Check the running of the cable.					
Check all remaining cables for damages. Also check the running of the cables.					
Check of siderails' fixing and locking/movement.					
Check the wheels (Security, braking and free rolling)					
The following parts must be lubricated: Centres of rotation (Motors and mattress support parts) Centres of rotation of rastofix fitting					

A service contract can be made in the countries, where Invacare[®] has its own sales company.

Moreover, Invacare[®] offers courses in service and maintenance of the Etude bed.

14. Waste disposal

All wooden parts must be dismantled and sent for incineration.

All electric parts must be dismantled and be disposed of as electric components.

Plastic parts must be sent for incineration or recycling.


Steel parts and castors must be disposed of as waste metals.

15. Trouble-shooting with the electrical system


Service and maintenance of the Etude must only be performed by personnel who have received the necessary instruction or training.

Symptom	Possible cause	Remedy
Mains indicator does not light up	1) Mains not connected 2) Fuse blown 3) Control damaged	1) Connect mains 2) Change the fuse, if control is designed for external fuse change or replace the control 3) Replace the control
Mains indicator lights up, but motor does not run. The relays in the control give an audible click	1) Motor plug not pushed correctly into the control 2) Motor damaged 3) Control damaged	1) Push the motor plug correctly into the control 2) Replace the motor 3) Replace the control
Mains indicator lights up, but motor does not run. No relay sound is heard from the control	1) Control damaged 2) Control box damaged	1) Replace the control 2) Replace the control box
Control is in order except one direction on one channel	1) Control box damaged 2) Control damaged	1) Replace the control box 2) Replace the control
No motor noise or piston rod movement	1) Motor plug not pushed correctly into the control 2) Fuse in the control blown	1) Push the motor plug correctly into the control 2) Change the fuse, if control is designed for external fuse change or replace the control
The motor is running, but the spindle does not move	Gear wheel or spindle damaged	Replace the motor
The motor cannot lift full load	Motor damaged	
Motor noise, but no movement of piston rod		
Piston rod will only operate inwards and not back outwards	Safety nut has become activated	


16. Wiring


1. Assemble the bed and connect 230 V mains plug to the control box.


2. Connect the bed end motor cables (head and foot end) to the control box. The cables must hang loose. Cable number: 02770-150-1250.


3. Raise the bed to the top position.

4. Place the foot end motor cable in the hook by the foot end. The backrest must be in a horizontal position.

5. 4-sec. bed: The plug of the thigh motor must be directed through the opening of the supporter for the backrest motor and then connected to the control box. Cable number: 0277011-1250.

6. Place the head end motor cables in the four hooks by the head end. The backrest must be raised.

Exchange of control box and cables

ATTENTION: The control box must always be placed with the plugs face up.

1. The control box can be removed by removing the fastener clips from the backrest motor.

Fastener clips

2. Put the control box in a sideways position from the motor.

3. The control box and/or cables can now be engaged.

17. Technical specifications

All measurements are stated in cm. All angles are stated in degrees. All measurements and angles are stated without tolerances.

Invacare® reserves the right to change the stated measurements and angles.
(Illustration 11)

Voltage supply: 230 V ~ ±10%, 50 Hz.

Max. current input: 1 A.

Voltage output: 24 V ~ max. 70 VA

Intermittent (periodic motor operation): 10%, max. 6 min/h.

Protection class: IP 54.

Insulation class: II, type B.


The patient is not separated from the ground and the chassis.

Double insulated.

The bed is not provided with a mains switch, so the mains plug is the only separation from the mains.

= Basic: 180 kg } - Max. load (SWL)
Classic: 220 kg } (Patient + mattress + siderail + lifting pole + other equipment)

Sound level: 45-50 dB (A)

18. Weight

Basic bed end - 1 piece	17 kg
Classic bed end - 1 piece	24 kg
Mattress support, head end	22 kg
Mattress support, foot end, 2-sectioned	17 kg
Mattress support, foot end, 4-sectioned	22 kg
Siderail - 1 piece	7 kg
Lifting pole	7 kg
Mattress support extension (10 cm)	3,2 kg
Mattress support extension (20 cm)	5,5 kg

2-sectioned bed, complete, excluding accessories with Basic bed ends	73 kg
4-sectioned bed, complete, excluding accessories with Basic bed ends	78 kg

2-sectioned bed, complete, excluding accessories with Classic bed ends	87 kg
4-sectioned bed, complete, excluding accessories with Classic bed ends	92 kg

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INVACARE Deutschland GmbH
INVACARE LTD
EC-Høng A/S**