

## User Instruction for the Bariatric Toilet Surround X154



### 1. Introduction

This Toilet Surround supports clients to give them confidence as they sit and rise – it offers support and lessening of strain. The height is adjustable to cater for varying patients requirements. The padded arms improve the users grip. The rubber feet add extra grip. This product is suitable for intensive cleaning.

Product Dimensions	
Width between arms	546mm
Frame depth	470mm
Adjustable height (from-to)	654mm-807mm
Footprint	495mm x 635mm

### 2. Contents of the Box

- Toilet Surround (X154)
- Extending Legs & Pin Clips
- User Instructions (this leaflet)


### 3. Assembly Instructions

- Carefully unpack and examine the product to ensure it is complete.
- Before use check for any signs of damage e.g. cracked weld, flaking of paint, bent tube. If any signs occur please contact your supplier to exchange and refrain from using the product. **DO NOT ATTEMPT TO USE PRODUCT IF DAMAGED.**

### 4. Correct Use

#### Adjustable Height:

- Remove the pin clips from all leg uprights.
- Select the desired height of the frame by gently sliding out each extension leg from the frame (NOTE – this is a telescopic action).



## **User Instruction for the Bariatric Toilet Surround**

### **X154**

- When desired height is obtained, refit the pin clip into all of the legs making sure pin clip goes through each leg and its extension leg. The frame is now ready for use.

#### **5. Cleaning**

- Use mild detergent with warm water to wipe down framework on at least a weekly basis more frequently if required.
- Do not use an abrasive cleaner as this will damage the frame and nullify your warranty.

#### **Use of power/jet washing & industrial washing machines:**

- If either of the above is used to clean equipment care must be taken since high-pressure water can be penetrative.
- Care should be taken that frames are not placed too close to high-pressure water jets.
- After cleaning ensure any residual water is removed from the framework as it could cause permanent rusting of the frame (leave frames open in a warm environment for several hours to ensure all internal moisture evaporates).
- We cannot offer advice on industrial washing machines with regards to specific machine settings, as these will be dependant on the manufacturer's instructions.

#### **6. Maintenance Requirements**

- The product should be checked on a regular basis for damage or wear and tear – this includes flaking or cracking paint. Any undue signs of damage should be noted and reported to the toilet surround supplier in the first instance. Should damage be evident, then the user should refrain from using the product until advice has been sought from the supplier.

#### **7. Warning**

- This product is intended for users of a maximum of 250kg/40stone in weight. TO BE REVISED.
- We recommend a full risk assessment should be carried out before use.

#### **8. Warranty**

This product carries a 12 month warranty from the date of purchase against manufacturing defects. Misuse or damage is not covered. This warranty does not affect your statutory rights. As manufacturers we cannot specify expected product life as this is based on the level of usage.

#### **9. Disposal Instructions**

- Appropriate hygiene standards should be applied before disposal to reduce the risk of infection.
- Materials used are essentially non hazardous and can be disposed of appropriately.

#### **10. Contact Details**

We have attempted to cover most eventualities, however if you have any specific queries, please contact:

Cefndy Healthcare  
Cefndy Road  
Rhyl  
Denbighshire  
LL18 2HG  
Tel: 0044 1745 343877 (Quality Assurance Dept)  
E-mail: [cefndyengineering@denbighshire.gov.uk](mailto:cefndyengineering@denbighshire.gov.uk)